


# Direct RNA Sequencing: Technique and Applications


Joseph Won

# Traditional RNA Sequencing



- In classic RNA-seq, complementary DNA (cDNA) is required
- Requires fragments with a poly-A tail attached
- Therefore, portions of the of the transcriptome are not sequenced

# Transcriptomics using RNA-seq




To differentiate from rRNA, all portions without a poly(a) are ignored

PCR amplification can introduce significant biases

# Direct RNA Sequencing

- Nanopore Sequencing: reads individual DNA or RNA strands by passing them through tiny protein nanopores, detecting changes in electrical current as different nucleotides pass through
- Benefits for RNA:
  - No conversion to cDNA
  - Full length RNA reads
  - Detection of RNA modifications (i.e. methylation)



Library preparation for MinION sequencing

# Requirements and Costs:



- Equipment:
  - Sequencing device (MinION)
  - Library Prep Kit
- Software:
  - Minimap2 (alignment tool)
  - Nanopolish (variant calling tool)

## Nanopolish

nanopolish passing

Software package for signal-level analysis of Oxford Nanopore sequencing data. Nanopolish can calculate an improved consensus sequence for a draft genome assembly, detect base modifications, call SNPs and indels with respect to a reference genome and more (see Nanopolish modules, below).

## Users' Guide

Minimap2 is a versatile sequence alignment program that aligns DNA or mRNA sequences against a large reference database. Typical use cases include: (1) mapping PacBio or Oxford Nanopore genomic reads to the human genome; (2) finding overlaps between long reads with error rate up to ~15%; (3) splice-aware alignment of PacBio Iso-Seq or Nanopore cDNA or Direct RNA reads against a reference genome; (4) aligning Illumina single- or paired-end reads; (5) assembly-to-assembly alignment; (6) full-genome alignment between two closely related species with divergence below ~15%.



The MinION Flow Cell can generate up to 50 Gb of data for sequencing DNA, cDNA or native RNA in real-time.

\$1,999.00

## Direct RNA Sequencing Kit

SQK-RNA004



The PromethION Flow Cell can generate up to 290 Gb for sequencing DNA, cDNA or native RNA in real-time.

\$436,000.00



A sequencing kit optimised for sequencing native RNA with improved output and accuracy.

This is an Early Access product.

Product lead time: 1 week

US\$599.00

- 1 +

Add to basket

2 Early Access

# Applications

Article | [Open access](#) | Published: 19 December 2022

## **Nano3P-seq: transcriptome-wide analysis of gene expression and tail dynamics using end-capture nanopore cDNA sequencing**

[Oguzhan Begik](#), [Gregor Diensthuber](#), [Huanle Liu](#), [Anna Delgado-Tejedor](#), [Cassandra Kontur](#), [Adnan Muhammad Niazi](#), [Eivind Valen](#), [Antonio J. Giraldez](#), [Jean-Denis Beaudoin](#), [John S. Mattick](#) & [Eva Maria Novoa](#) 

*Nature Methods* **20**, 75–85 (2023) | [Cite this article](#)

19k Accesses | 18 Citations | 54 Altmetric | [Metrics](#)

Article | [Open access](#) | Published: 14 February 2019



## **Direct RNA sequencing on nanopore arrays redefines the transcriptional complexity of a viral pathogen**

[Daniel P. Depledge](#) , [Kalanghad Puthankalam Srinivas](#), [Tomohiko Sadaoka](#), [Devin Bready](#), [Yasuko Mori](#), [Dimitris G. Placantonakis](#), [Ian Mohr](#) & [Angus C. Wilson](#) 

*Nature Communications* **10**, Article number: 754 (2019) | [Cite this article](#)

32k Accesses | 138 Citations | 60 Altmetric | [Metrics](#)

## **Identification of differential RNA modifications from nanopore direct RNA sequencing with xPore**

[Ploy N. Pratanwanich](#) , [Fei Yao](#), [Ying Chen](#), [Casslynn W. Q. Koh](#), [Yuk Kei Wan](#), [Christopher Hendra](#), [Polly Poon](#), [Yeek Teck Goh](#), [Phoebe M. L. Yap](#), [Jing Yuan Chooi](#), [Wee Joo Chng](#), [Sarah B. Ng](#), [Alexandre Thiery](#), [W. S. Sho Goh](#)  & [Jonathan Göke](#) 

*Nature Biotechnology* **39**, 1394–1402 (2021) | [Cite this article](#)

23k Accesses | 128 Citations | 134 Altmetric | [Metrics](#)

- **Gene expression analysis:**  
Capturing full-length transcripts including exons, introns, and untranslated regions
- **Splicing and isoform diversity:**  
Detecting alternative splicing events
- **Epitranscriptomics:**  
Identification of RNA modification directly on native RNA molecules



# Limitations


**Error rates:** Higher error rates compared to DNA sequencing

“Since 2019, two publications have reported median accuracy of ~91% for *Brassica napus* and ~88–90% for *E. coli*”

**Lower throughput:** Therefore, lower depth in sequencing and requires higher input RNA quantity

**Less computational tools:** By virtue of its novelty, there are less tools and datasets to utilize nanopore sequencing data

# Question: Where do you see application of dRNA sequencing in evolutionary biology?



Keep in mind MLS is controversial because there is no good way to link different “levels” of biological organization – even though we can obviously observe the patterns (differential gene expression, phenotypes, social behavior, etc...)

## Pattern vs. Process

### Traditional RNA Seq:

- Gene expression quantification
- Differential gene expression

### Direct RNA Seq:

- Isoform Diversity
- RNA Modifications
- Methylation
- Alternative Splicing

\*as related to the transcriptome